


## Features

The number of output channels of ROI - SCA is expanded to 4, and TTL (pulse width 50 ns) can be output even at FAST timing (setting example 1). Also, you can assign ROI - SCA functions to the expanded 4 channels (setting example 2).



Setting example 1

AUX	type of AUX signal
AUX1	ROI1-SCA
AUX2	ROI2-SCA
AUX3	ROI3-SCA
AUX4	ROI4-SCA
AUX5	fast-CH1
AUX6	fast-CH2
AUX7	fast-CH3
AUX8	fast-CH4

Setting example 2

AUX	type of AUX signal
AUX1	ROI1-SCA
AUX2	ROI2-SCA
AUX3	ROI3-SCA
AUX4	ROI4-SCA
AUX5	ROI5-SCA
AUX6	ROI6-SCA
AUX7	ROI7-SCA
AUX8	ROI8-SCA

## Expand count upper limit

Extend the count number of conventional quick scan from 16 bit to 32 bit, and support step scan

## Features

We have expanded the upper limit of the number of quick scan counts (high-speed histogram acquisition function) supported by the APN (U) 504 series. It is possible to correspond to step scan which can be measured for a long time. Quick scan data can be expanded by adding input counts.

### Conventional Quick scan data structure (1 event: 32768byte)

Event#1	CH1 4096ch 16bit	CH2 4096ch 16bit	CH3 4096ch 16bit	CH4 4096ch 16bit
Event#2	CH1 4096ch 16bit	CH2 4096ch 16bit	CH3 4096ch 16bit	CH4 4096ch 16bit

### Extended histogram depth from 16 to 32-bit (1 event: 65536 byte)

Event#1	CH1 4096ch 32bit	CH2 4096ch 32bit	CH3 4096ch 32bit	CH4 4096ch 32bit
Event#2	CH1 4096ch 32bit	CH2 4096ch 32bit	CH3 4096ch 32bit	CH4 4096ch 32bit

### Extended by adding input count data (1 event: 65552 byte)

Event#1	CH1 4096ch 32bit	CH2 4096ch 32bit	CH3 4096ch 32bit	CH4 4096ch 32bit	CH1 Input 32bit	CH2 Input 32bit	CH3 Input 32bit	CH4 Input 32bit
Event#1	CH1 4096ch 32bit	CH2 4096ch 32bit	CH3 4096ch 32bit	CH4 4096ch 32bit	CH1 Input 32bit	CH2 Input 32bit	CH3 Input 32bit	CH4 Input 32bit

\*Images is for illustration purpose.  
\*Please note that contents may change without prior notice.

Website



Manufacture of Radiation and Radioactivity measurement devices

# TechnoAP Co., Ltd.

2976-15 Mawatari, Hitachinaka-shi, Ibaraki, 312-0012, Japan

+81-29-350-8011

+81-29-352-9013

order@techno-ap.com