

APG7300D USB-MCA

Command Manual

Version 1.0

May ,2025

TechnoAP Co., Ltd.

ADD : 2976-15 Mawatari, Hitachinaka-shi, Ibaraki, Japan

ZIP Code: 312-0012

TEL: 029-350-8011

FAX: 029-352-9013

URL: <http://www.techno-ap.com>

e-mail: order@techno-ap.com

Content

1 Safety Precautions / Disclaimer	3
2 Overview	4
2.1 Overview	4
2.2 Revision History	4
2.3 Environment.....	4
3 Commands	5
3.1 Overview	5
3.2 Command List.....	8
3.3 Command Description	9

1 Safety Precautions / Disclaimer

Thank you very much for purchasing this product from TechnoAP Co., Ltd. Before using this product, please read this "Safety Precautions and Disclaimer" and be sure to observe the contents and use the product properly.

We are not responsible for any damage caused by abnormality of device, detector, connected device, application, damage to failure, other secondary damage, even if accident caused by using this device.



Prohibited matter

This device cannot be used for applications requiring special quality and reliability related to human life, accident.

This device cannot be used in places with high temperature, high humidity, and high vibration.

Do not apply a power supply that exceeds the rating.

Do not turn the power on while other metals are in contact with the board surface.



Note

If there is smoking or abnormal heat generation in this device, turn off the power immediately.

This board may not work properly in noisy environments.

Be careful with static electricity.

The specifications of this board and the contents of the related documents are subject to change without notice.

Warranty policy

The warranty conditions of "our product" are as follows.

Warranty period: One year from date of purchase.

Guarantee contents: Repair or replacement will be carried out in case of breakdown even though you have used correctly according to this instruction manual within the warranty period.

Out of warranty: We do not warranty if the cause of the failure falls under any of the following.

1. Failure or damage due to misuse or improper repair or modification or disassembly.
2. Failure and damage due to falling etc.
3. Consumables.

2 Overview

2.1 Overview

This device can be connected to a PC via a USB cable and can perform measurement control by sending and receiving commands from environments such as Microsoft Visual Studio, National Instruments LabVIEW, and Linux.

This document explains the handling of measurement control commands for this device.

The contents described in this document are subject to change without prior notice.

2.2 Revision History

May 12, 2025 Version 1.0 First edition

2.3 Environment

The following environment is assumed to be prepared.

- (1) FTDI D2XX Drivers are installed.
- (2) Our sample program can be executed.

3 Commands

3.1 Overview

All settings and data acquisition for this device are performed via USB communication.

The commands are broadly classified into **Setting**, **Status**, **Histogram**, and **List** commands. The details of each type are described below.

(1) Setting Commands

For setting commands, the PC sends 8 bytes, and the device returns an 8-byte response.

After sending the 8 bytes, the PC must immediately read the 8-byte response.

Example: Setting a mode

PC →Device	Command field (4 bytes) ASCII string 4D4F4457 (MODW)	Parameter field (4 bytes) Binary 0x00000000		
			Command field (4 bytes) ASCII string 4D4F4457 (MODW)	Parameter field (4 bytes) Binary 0x00000000
Device →PC				

[Configuration from the PC]

The command field is 4 bytes and is an ASCII string.

For example, when setting the mode, it is “MODW”, which corresponds to “4D4F4457” in ASCII representation.

The parameter field is a 4-byte binary value in big-endian format (network byte order, MSB first) and represents the value to be set.

[Response from the device]

If the configuration is executed correctly, a response with the same contents as the configured value is returned.

By comparing the configured value with the response value, it is possible to confirm whether the configuration command was executed correctly.

(2) Status Commands

For status commands, when the PC sends 4 bytes, the device returns a 25-byte response.

After sending 4 bytes, it is necessary to immediately read the 25-byte response.

PC →Device	Command field (4 bytes) ASCII string 53545557(STUW)	Parameter field (4 bytes) Binary 0x00000000		
Device →PC			Data field (25Byte) Binary real time, CH1 live time, dead time, throughput rate, throughput count	

[Setting from the PC]

The command field is 4 bytes and is an ASCII string.

It is "STUW", which corresponds to "53545557" in ASCII representation.

The parameter field is a 4-byte binary value and is set to 0.

[Response from the device]

If the command is executed correctly, real time, CH1 live time, dead time, throughput rate, and throughput count are received together.

(3) Histogram Commands

Histogram data has a maximum size of 65536 bytes per channel(65536 bytes = 16384 channels × 4 bytes per channel). Since reading is performed in units of one block (2048 bytes), it is necessary to divide the read operation while specifying the block number to be read.

When the PC sends a histogram read command including block number information (8 bytes), the device returns a 2048-byte response.

After sending the read command, it is necessary to immediately read the 2048-byte response.

This 2048-byte response does not include a header and contains only the histogram data for the specified block.

Each block contains data for 512 channels


(512 = 2048 bytes / 4 bytes).

To acquire the maximum of 16384 channels, it is necessary to execute the operation 32 times


(32 = 16384 channels / 512 channels).

Example: Acquiring histogram data for CH1 for 16384 channels

Histogram data read



* Continue executing in the same manner with block numbers 1 through 31.



[Histogram read and block number setting from the PC]

The command field is 4 bytes and is set as an ASCII string in the "HIXX" format, where XX represents the block number to be read in hexadecimal.

To acquire channels 0 to 511, the ASCII string is "HI00".

To acquire channels 15872 to 16383, the ASCII string is "HI1F".

The parameter field is a 4-byte binary value and is set to 0.

The histogram data for the specified channel range of the selected channel is sent from the device as 2048 bytes corresponding to 512 channels.

The same procedure is repeated while specifying ranges as needed to read the required number of channels.

3.2 Command List

Command			Description	Remarks		
Category	Abbreviation	ASCII				
Register Setting	CH1(ch.0)	ADGW	41444757	CH1 ADC Gain		
		THRW	54485257	CH1 Threshold		
		LLDW	4C4C4457	CH1LLD		
		ULDW	554C4457	CH1ULD		
		OFSW	4F465357	CH1 Offset		
	Common	MODW	4D4F4457	Mode	0: Histogram (fixed)	
		MMDW	4D4D4457	Measurement Mode	0 : real time, 1 : live time	
		MT0W	4D543057	Measurement Time	Upper 12 bits	
		MT1W	4D543157		Lower 32 bits	
		PDSW	50445357	Peak Detection Mode		
		AQSW	41515357	Measurement Start		
		AQEW	41514557	Measurement Stop		
		CLRW	434C5257	Clear		
	Status		STUW	53545557	Status	real time, CH1 live time, etc.
	Histogram	HI00	48493030	Histogram Memory Read Read in units of 512 points × 4 bytes	0~511Channels	
HI01		48493031	512~1023 Channels			
HI02		48493032	1024~1535 Channels			
HI03		48493033	1536~2047 Channels			
HI04		48493034	2048~2559 Channels			
HI05		48493035	2560~3071 Channels			
HI06		48493036	3072~3583 Channels			
HI07		48493037	3584~4095 Channels			
HI08		48493038	4096~4607 Channels			
HI09		48493039	4608~5119 Channels			
HI0A		48493041	5120~5631 Channels			
HI0B		48493042	5632~6143 Channels			
HI0C		48493043	6144~6655 Channels			
HI0D		48493044	6656~7167 Channels			
HI0E		48493045	7168~7679 Channels			
HI0F		48493046	7680~8191Channels			
HI10		48493130	8192~8703 Channels			
HI11		48493131	8704~9215 Channels			
HI12		48493132	9216~9727 Channels			
HI13		48493133	9728~10239Channels			
HI14		48493134	10240~10751Channels			
HI15		48493135	10752~11263Channels			
HI16		48493136	11264~11775Channels			
HI17		48493137	11776~12287Channels			
HI18		48493138	12288~12799Channels			
HI19		48493139	12800~13311Channels			
HI1A		48493141	13312~13823Channels			
HI1B		48493142	13824~14335Channels			
HI1C		48493143	14336~14847Channels			
HI1D		48493144	14848~15359Channels			
HI1E		48493145	15360~15871Channels			
HI1F		48493146	15872~16383Channels			

3.3 Command Description

*For details of each setting, refer to the “USB-MCA User’s Manual” supplied with the device.

(1) ADC Gain

Description	Sets the ADC gain and the number of channels.
Range	0: 16384 channels 1: 8192 channels 2: 4096 channels 3: 2048 channels 4: 1024 channels 5: 512 channels

(2) Threshold

Description	Sets the threshold level for waveform acquisition. The value must be equal to or lower than the LLD setting.
Range	0 to 16383

(3) LLD

Description	Energy LLD (Lower Level Discriminator). The unit is channel count. The value must be equal to or greater than the threshold and less than the ULD.
Range	0 to 16383

(4) ULD

Description	Energy ULD (Upper Level Discriminator). The unit is channel count. The value must be greater than the LLD.
Range	0 to 16383

(5) Offset

Description	Sets the positive offset value in channel units.
Range	0 to 16383

(6) Mode

Description	Sets the operation mode.
Range	0: Histogram (fixed)

(7) Measurement Time

Description Sets the measurement time.

The internal clock operates at 25 MHz, so the value is calculated as:

Measurement time (seconds) × 25,000,000r

The value is divided into two registers: upper 12 bits and lower 32 bits.

Range Maximum 192 hours

(60 seconds × 60 minutes × 192 hours × 25,000,000 = 17,280,000,000,000)

In this case, set 0xFB7 to the upper register and 0x50430000 to the lower register.

(8) Peak Detection Method

Description Selects the method for detecting the peak (maximum pulse height).

Range 0: abs

1: fast

(9) Measurement Start

Description Starts the measurement.

Range 1: Start

(10) Measurement Stop

Description Stops the measurement.

Range 1: Stop

(11) Clear

Description Clears the measurement time and measurement data.

Range 0: Clear

(12) List Data Read Start

Description Starts reading list data.

Range 0 (fixed)

(13) Status

Description Retrieves status information.

No.	Item Description	Size (Byte)
1	Real Time (40 ns per count)	6
2	CH1 Live Time (40 ns per count)	6
3	CH1 Dead Time (40 ns per count)	6
4	CH1 Throughput Count Rate	3
5	CH1 Throughput Total Count	4
Total		25