



We specialize in radiation measuring devices

Allow data collection at extremely high rates with stable peak position and energy resolution



High throughput data acquisition


Readout for low dead time

APV8104-14

- **4CH 1GHz, 14bit**
- **Coarse: 2ns | Fine: 7.8ps**
- **Throughput 1Mcps and more / CH**
- **Functions: Digital CFD, TDC and QDC**

Please visit below for more information.

<http://www.techno-ap.com>



CUSTOM MADE for your needs. Contact us!

TechnoAP Co., Ltd.

2976-15 Mawatari, Hitachinaka-shi, Ibaraki, 312-0012

Manufacturer of Radiation Measurements

☎ +81-29-350-8011

☎ +81-29-352-9013